EN: Instructions for use

Only for use with the Sepax™ C-Pro and Sefia™ S-2000.

Before use:

Check pouch/blister is sealed and sterile. Check expiry date. Check for twisted, crushed, or kinked tubing. Do not use if defective. Inform Cytiva and return the kit.

Important:

This accessory must be used with a sterile connection device (SCD).

Instructions for use

1. Pouch can be opened outside of aseptic conditions.
2. Close the clamp on the sampling line, approximately 2 cm (1") from the sampling chamber.
3. Connect to the input bag with an SCD. Leave 20 cm of tubing to the input bag.
4. Verify the integrity of the connection.
5. Position the input bag vertically.
6. Open the clamp on the sampling line.
7. Compress and release the sample chamber, to draw product. Repeat until required volume is reached.
8. For the CPAK-100 only: if applicable, draw the product using Sample Calibrator accessory.
9. Invert the position of the input bag.
10. Compress and release the sample chamber, to push the residual product in the line back into the input bag.
11. Close the clamp on the sampling line.
12. Using a validated procedure: Seal the line near the sampling chamber. Detach the sampling chamber. Transfer the sample to a syringe under laminar flow.

Note: The volumetric scale on the sampling chamber provides an indicative volume only.

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cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Sepax and Sefia are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2020–2024 Cytiva

For local office contact information, visit cytiva.com/contact