## EN: Instructions for use

**Only for use with the Sepax™ C-Pro.**  
**Before use:**

Check pouch/blister is sealed and sterile. Check expiry date. Check for twisted, crushed, or kinked tubing. Do not use if defective. Inform Cytiva and return the kit.

### Instructions for use

- Under aseptic conditions, open the blister and close the input line clamp.
- Connect the initial product bag to the input line with an SCD outside laminar flow. Connect under aseptic conditions if you use the spike.
- Connect output bag to output line with SCD outside laminar flow. Connect under aseptic conditions if you use the Luer connector or the spike.
- Check the stopcock position is TTT.
- Install the kit on Sepax C-Pro system, following the instructions in the *Operator Manual*.

This kit is compatible with 4-6 mm membrane ports.

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**cytiva.com**

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Sepax is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2020–2024 Cytiva

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)