EN: Instructions for use

Only for use with the Sepax™ C-Pro.
Before use:

Ensure the pouch/blister is sealed and sterile. Check the expiry date. Check for twisted, crushed, or kinked tubing. If defective, do not use. Inform Cytiva and return the kit.

Important:

Perform all connections under aseptic conditions or use a sterile connection device (SCD).

Instructions for use

- Under aseptic conditions, open the blister and close all clamps, including those on the extension lines.
- Connect the washing solution bag to the spike at top right. (Red stopcock.)
- Connect the output bags (one or two output bags, depending on requirements) to the spikes at lower left. (Blue stopcock). Seal any unused line.
- Check the stopcock position is TTT.
- Install the kit on Sepax C-Pro system, following the instructions in the *Operator Manual*.
- Connect the input bag using the appropriate connection extension with SCD. Seal any unused line.

This kit is compatible with 4-6 mm diameter ports.

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cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Sepax is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2020–2024 Cytiva

For local office contact information, visit cytiva.com/contact