## EN: Instructions for use

**Only for use with the Sepax™ C-Pro.**  
**Before use:**

Ensure the pouch/blister is sealed and sterile. Check the expiry date. Check for twisted, crushed, or kinked tubing. If defective, do not use. Inform Cytiva and return the kit.

### Important:

Perform all connections under aseptic conditions or use a sterile connection device (SCD).

### Instructions for use

- Under aseptic conditions, open the blister and close all clamps, except the red one.
- Connect a 20 mL syringe to the Luer port. If not required, seal the unused line.
- Connect the washing solution bag to the spike at top right. (Red stopcock).
- Connect the re-suspension medium bag to the spike at top right. If not required, seal the unused line.
- If required, fill the bag at bottom right with the required amount of density gradient medium. Use the filtered port. Seal near the filtered port when done.
- Connect the initial product bag to the input line.
- Check the stopcock position is TTT.
- Install the kit on Sepax C-Pro system, following the instructions in the *Operator Manual*.

This kit is compatible with 4-6 mm diameter ports.

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**cytiva.com**

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Sepax is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2020–2024 Cytiva

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)