EN: Instructions for use

Only for use with the Sefia™ platform.

Before use:

Ensure the pouch/blister is sealed and sterile. Check the expiry date. Check for twisted, crushed, or kinked tubing. If defective, do not use. Inform Cytiva and return the kit.

Important:

Perform all connections under aseptic conditions or use a sterile connection device (SCD).

Instructions for use

1. Open the blister.
2. Close all clamps.
3. If outside aseptic conditions, install the kit on the machine and use an SCD to connect the initial product to the initial line. If under aseptic conditions, connect the initial product via the spike and then install the kit on the machine. See the *operator manual* for installation instructions.
4. Check the stopcocks are correctly aligned. See the kit section in the *operator manual* for the correct position.

This kit is compatible with 4-6 mm diameter ports.

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cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Sefia is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2020–2024 Cytiva

For local office contact information, visit cytiva.com/contact