

Dispenser pack refill for Sentino™ filter dispenser

Product information

Introduction

Important

Read these instructions carefully before using the products.

Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

Description

This dispenser pack refill is designed for use with the Sentino™ filter dispenser.

When loaded into the automatic filter dispenser, individual sterile membrane filters can be aseptically dispensed. During dispensing, the packaging waste winds onto spools and can be discarded at the end of the pack.

Specifications

The specifications for each membrane type are listed below.

Membrane	Supor™	ME25	Metrcel™ Black	
Pore size	0.2 μm	0.45 μm	0.45 μm	0.8 μm
Media type	Polyethersul-fone (PES)	Mixed Cellulose Ester (MCE)	Modified PES	
Recovery compared to spread plates	Meets or exceeds global recovery performance standards. Note: <i>Recovery rates and colony morphology of test organisms indicate there is no enhancement or inhibition due to:</i> • Chemical extractables. • Grid lines. • Gamma irradiation. • Packaging materials.			
Test organism	<i>P. aeruginosa</i>	<i>E. coli</i>	<i>S. cerevisiae</i>	
Sterilization	Gamma-irradiated			

Instructions for use

For loading instructions refer to *Sentino filter dispenser Operating Instructions (89965)*.

Ordering information

For additional ordering information, visit [cytiva.com](https://www.cytiva.com).

Product code	Description	Pkg ^{1,2}
68121ME	ME25 MCE membrane 0.45 µm	1000/pkg
68123	Supor PES membrane 0.2 µm	1000/pkg
68124	Metrcel Black PES membrane 0.45 µm	1000/pkg
68125	Metrcel Black PES membrane 0.8 µm	1000/pkg

¹ Each box contains two additional discs for loading the dispenser.

² Each pack of 1000 contains 5 boxes, each box contains 200 filters.



Give feedback on this document

Visit cytiva.com/techdocfeedback or scan the QR code.



cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Metrical, Sentino, and Supor are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Any other third-party trademarks are the property of their respective owners.

© 2024–2025 Cytiva

For local office contact information, visit cytiva.com/contact

29743967 AC V:7 06/2025

